



FURIX

FVP-60



Furix Veterinary Infusion Pump



13 Infusion Modes

Rate Mode, Rate Time Mode, Rate Volume Mode, Body Weight Mode, Drug Library Mode, Sequential Mode, Loading dose Mode, Micro Mode, Ramp Up/Down Mode, Intermittent mode, Relay Mode, Drip Mode, TIVA Mode(Optional)

5+

Open/Close system Option

Open system match different brands of infusion set with self calibration to guarantee best infusion accuracy
Close system match special infusion set to guarantee best infusion accuracy



UI Design

Friendly and clear for monitoring of infusion status with Rate, time-left, VTB, Accumulated volume, Pressure, Drug name etc



DPS Function

DPS(Dynamic pressure detection) function, monitoring the pressure variations in real time



Big Touch Screen

4 inch colorfu & touch screen for easy operation



User Customization

Drug library & Syringe brand can be modified by the touchscreen keyboard



Upper/Downforce Occlusion Detection

Make the injection more secure



System Self Test

Safeguard for every time use

Technical Specification

Function	Mode	FVP-60
LCD		4 Inch Touch Screen
Infusion Mode		Rate Mode, Rate Time Mode, Rate Volume Mode, Body Weight Mode, Drug Library Mode, Sequential Mode, Loading dose Mode, Micro Mode, Ramp Up/Down Mode, Intermittent Mode, Relay Mode, Drip Mode, TIVA Mode(Optional)
Infusion Principle		Fingertip Peristaltic Pump
I.V Set Compatibility		Open System, match all brands of qualified PVC, TPE I.V sets with external diameter of 3.8mm 2mm
Range of Flow Rate		0.01-2200ml/h
Flow Rate Step by		0.01ml/h
Purge/Bolus Rate		1-2200ml/h, default 800ml/h, step by 0.01ml/h
Single Bolus Volume		0.1-50ml Adjustable, default 3ml
Time Preset		00:01~99:59 (Hour: Min)
Volume Range		1~9999ml
Drop Range		1-600drops/min, step by 1 drop
Accuracy		±5%
Airline Detection		Ultrasonic, Upper and under detection
Visible & Audible Alarms		Air in line, Occlusion, Door Open, VTB Completion, Near Completion, Low Battery, Battery Exhausted, Low Temperature, Motor Malfunction, A/C Disconnection, A/C connection, I.V set dislocation, MPU Error, Battery Charging, Battery Charging Completion, Circuit Malfunction
KVO		1ml/h 5ml/h, default value 1ml/h, can be programmed by user step 0.1ml/h
Occlusion Pressure		Low 0KPa ± 20KPa Middle 60KPa ± 20KPa High 100KPa ± 20KPa Can be customized to 8-13 levels
Internal Battery		Lithium battery, 10.8/2200mAh, more than 8 hours working back up
Power Supply		30VA; AC100V 2.0V 50HZ/60HZ
Classification		Class II, type CF, PX
Dimension & Weight		22.16 x 6.9cm, 2KG
Optional Function		Ambulance DC12V, W/F, Language customization

Shenzhen Falcon and Furix Lifesciences Co. Ltd
Furix Trade Venture Limited (HK)

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
Mob: +852 9579 0472 E-mail: sales@myfurix.com

FURIX