



FURIX

FVM12 PRO

Multiparameter Patient Monitor



Product specification

	ECG	SPO2	(NIBP)	(RESP)	(TEMP)	HR/PR	(ETCO2)	(IBP)
Standard	✓	✓	✓	✓	✓	✓		
Optional							✓	✓

Function Configuration :Central workstation/Available with trolley and wall bracket/Built-in printer

Product Features :

- 12-inch high-brightness LCD monitor.
- Mirror display, new visual effect, more intimate protection monitor.
- Multiple working modes: diagnostic, monitoring, surgical mode.
- Suitable for different groups of pets: large, medium, small pets.
- Multiple display modes: standard, large character, trend graph, blood pressure list, etc.
- Support multi-lead display of ECG waveforms on the same screen.
- Equipped with a blood pressure cuff for pets Accurate blood pressure measurement.
- Multiple analysis functions.
- Quick Setup Alarm Function.
- Built-in high-capacity lithium battery, long standby time.
- Suitable for use in a variety of conditions such as operating rooms.

Accessories for animals :



Blood pressure extension tubing



ECG leadwire



Oximetry Probes



Oxygen Extension Cord



Body temperature sensors



Animal Cuffs

Parameter monitoring applicable to animals

- NIBP: Anti-weak signal and anti-interference algorithms designed for the characteristics of animal blood flow, enabling fast and accurate blood pressure measurements.
- SpO2: advanced anti-interference and weak perfusion technology
- ECG: Patented ECG algorithm, supports multiple arrhythmia identification
- Mainstream and bypass CO2: ideal for measurement in spontaneously breathing and mechanically ventilated animals
- Efficient and easy to use

Multi-scenario test

Adapt to changing application environments



0.75 meter drop resistance



Resistant to a wide range of cleaning and sanitizing agents

Low-power design

Power saving and durable



Fan-less design



Low electricity consumption



Long life of the original



Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG
 E-mail: vet.furix@gmail.com, sales@myfurix.com Web: www.myfurix.com
 Tel: +852 9579 0472 Mob: +852 9579 0472

FURIX